



## TO WHOM IT MAY CONCERN

---

This letter confirms that Medisafe Sonic Irrigator PCF System has been tested and is validated to clean *da Vinci EndoWrist* instruments. For complete and current listing of *da Vinci EndoWrist* instruments validated with the SI PCF system, visit Intuitive Surgical website at [www.intuitivesurgical.com](http://www.intuitivesurgical.com).

The Medisafe SI PCF system meets ISO15883 / HTM2030, ASTM E2314-0, AAMI TIR30 and AAMI TIR12 standards. The SI PCF system has been tested to reprocess up to 10 *EndoWrist* instruments through its dedicated *da Vinci* cycle and uses the same critical parameters as outlined in the Intuitive Surgical's manual cleaning method.

When cleaning *da Vinci EndoWrist* instruments, it is important that the hospital follows the cleaning methods outlined in the Instrument Instructions For Use (IFU) and/or IFU Supplements that come in the instrument package.

When considering automated cleaning option for cleaning *da Vinci* instrument, it is important that the hospital uses only validated machines and requests a copy of the validation certificate from any vendor claiming the cleaning compatibility.

For any further information or details about the manual cleaning or automated cleaning options to clean *da Vinci EndoWrist* instruments, please contact Intuitive Surgical Customer Service at (888) 409-4774 or (408) 523-2100 or visit Intuitive Surgical website at [www.intuitivesurgical.com](http://www.intuitivesurgical.com).