

Certificate

SMP GmbH
Pruefen Validieren Forschen

hereby certifies the

cleaning validation of Ratched Instruments
(33 cm and 45 cm; rotatable and rigid)

of the company

Styker Endoscopy
5900 Optical Court
San Jose, CA 95138
USA

The SMP investigation, with the project no. 05108012604, shows that the cleaning effectiveness of the Medisafe SI PCF Washer in combination with the detergent 3 E Zyme is powerful enough to assure good cleaning results of Stryker Endoscopy Ratched Instruments without the need of manual pre-cleaning:

- All instruments are visibly clean
- All tips of the instruments are clean after the process.
- All tests fulfil the acceptance criteria

Tests followed

TS/ISO 15883-5
AAMI TIR 30
AAMI TIR 12



Tübingen, Germany; June 30. 2008

A handwritten signature in black ink, appearing to read 'K. Roth', is written over a horizontal line.

Klaus Roth
Managing Director